

There are no translations available.



M 45 – Pleiades

Date/time:

15.01.2009 18:55 to 15.01.2009 19:25

25.01.2009 17:38 to 25.01.2009 19:00

26.01.2009 18:00 to 26.01.2009 18:30

30.01.2009 18:00 to 30.01.2009 18:30

Location: Aach, Germany

Equipment:

Telescope: TMB 92 SS f/5.5 with TSFlat

Mount: CGE, selfguided with ST-4000XCM

Camera: SBIG ST-4000XCM

Filters: IDAS LPS

Acquisition:

Software: CCDSoft V5

FOV, Image scale: 99.4 x 94.5 arcmin, 3.01 arcs/pixel

Sub-exposures: 2h 30min (3 x 10min + 8 x 10min + 1 x 30min + 3 x 10min) light frames
@-35°C

45 x 10min, 37 x 30min dark frames

30 x 0,5s EL flats

Processing:

Calibration, stacking and Exc Avg combine of all sets in Images Plus 3.75,

curves and levels using adjustment layers in CS3,

star saturation increase with Carbonis Actions,

color enhancement using Gabany method under CS3,

green cast elimination with Pixinsight LE 1.0,

gradient reduction with GradientXterminator,

final noise reduction of background and sharpening of nebulosity in CS3

Comments:

These were the best imaging sessions on M45.

This image was a first attempt to stack frames from different sessions in IP 3.75.

{downloader mode="level" level="public" file="clusters/m45_2009-01-30_full.jpg"}